

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B067	TRICHOPHYTON AGAR-7				
Formula					
Ingredients :		gms/lit.			
Ammonium nitrate		1.50			
Dextrose		40.00			
Monopotassium dihydrogen phosphate		1.80			
Magnesium sulphate		0.10			
L-Histidine hydrochloride		0.03			
Nicotinic acid		0.002			
Agar		15.00			
Final pH (at 25°C) : 6.8 ± 0.2					
Directions :					
Suspend 58.43 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Dispense in test tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Allow the tubed medium to cool in a slanted position.					
Principle :					
The method also employs an ammonium nitrate basal medium (Trichophyton Agar-6, B068) to which histidine is added.					
QC Tests - (I) Dehydrated Medium					
	Colour :	White to light yellow			
	Appearance :	Homogeneous Free Flowing powder			
(II) Rehydrated medium					
	pH (post autoclaving/heating) :	6.8 ± 0.2			
	Colour (post autoclaving/heating) :	Light amber			
	Clarity (post autoclaving/heating) :	clear to slightly opalescent gel forms in tubes as slants			
(III) Q.C. Test Microbiological					
Cultural characteristics observed after an incubation at 25-30°C for 1 weeks.					
	MICROORGANISM (ATCC)	GROWTH			
	Microsporium gallinae (12108)	good-luxuriant			
	Trichophyton megninii (12106)	good-luxuriant			
Precautions :	1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :	1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :	It is used for differentiation of Trichophyton species				
Storage :	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :	500 gm. bottle				
Product profile:	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B067	58.43 g/l	8.557 L	6.8 ± 0.2	Nil	121°C /15 min.

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.