

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B1133		KANAMYCIN ESCULIN AZIDE AGAR BASE				
Formula						
Ingredients :		gms/lit.				
Tryptone		20.00				
Yeast extract		5.00				
Sodium chloride		5.00				
Sodium citrate		1.00				
Esculin		1.00				
Ferric ammonium citrate		0.50				
Sodium azide		0.15				
Agar		10.00				
Final pH (at 25°C) : 7.0 ± 0.2						
Directions :						
Suspend 21.32gms in 500ml distilled water. Add one vial of Kanamycin sulphate selective supplement(BF130) reconstituted as directed. Bring to boil. Dispense into final containers & Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.						
WARNING: Sodium azide has a tendency to form explosive metal azides with plumbing materials. It is advisable to use enough water to flush off the disposables.						
Principle :						
Tryptone, yeast extract provides essential nutrients for Enterococci. Kanamycin sulphate and sodium azide are the selective inhibitory components. Esculin and ferric ammonium citrate together form indicator system to detect esculin – hydrolysing Streptococci forming black zones around the colonies.						
QC Tests – (I)Dehydrated Medium						
	Colour :	Cream to light yellow				
	Appearance :	Homogeneous Free Flowing powder				
(II)Rehydrated medium						
	pH (post autoclaving/heating) :	7.0 ± 0.2				
	Colour (post autoclaving/heating) :	Medium amber				
	Clarity (post autoclaving/heating) :	Slightly opalescent gel with purpish tinge forms in petri plates.				
(III)Q.C. Test Microbiological						
	Cultural characteristics observed after	18-24 hrs at 35-37°C or 42°C.				
	MICROORGANISM (ATCC)	GROWTH	ESCULIN HYDROLYSIS			
	Enterococcus faecium (19434)	Luxuriant	+			
	Streptococcus bovis (27960)	Luxuriant	+			
	Escherichia coli (25922)	Inhibited	-			
Precautions :		1. For Laboratory Use.				
		2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		For selective isolation and identification of group D Streptococci in foodstuffs.				
Storage :		Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B1133		42.64 g/l	11.726L	7.0 ± 0.2	Kanamycin sulphate selective supplement (BF130)	121°C / 15 minutes

Refer disclaimer Overleaf

Rev: December 2020

BIOMARK Laboratories-INDIA
www.biomarklabs.com
TECHNICAL SHEET

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.