BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

B1455 YEAST EXTRACT AGAR						
Formula						
Ingredients:			gms/lit.			
Truntere			6.00			
Tryptone Yeast extract			3.00			
			15.00			
Agar	<u>``. 72 . 02</u>		12.00			
Final pH (at 25°C Directions :	.): 7.2 <u>+</u> 0.2					
	a in 1000 ml diati	llad wata	- llast to be	iling to discal	va tha madium	annalataly Ctariliza
Suspend 24 gram			r. Heat to be	ning to dissol		completely. Sterilize
	15 lbs pressure (121°C) fC	or 15 minute	15.001 to 45	-50°C. Mix well	and pour into sterile
Petri plates. Principle :						
	h nutrianta ara r	manidad	by twotopo	and vesst	avtraat Thaaa	
Necessary growth nutrients are provided by tryptone and yeast extract. These serve as source of						
nitrogen, vitamins, growth factors as well as crude source of carbon. Agar acts as a solidifying agent.						
2C Tests – (I)Dehydrated Medium			Our and the light wellow			
Colour :			Cream to light yellow			
			Homogeneous Free Flowing powder			
(II)Rehydrated medium			7.2.4.0.2			
			7.2 ± 0.2			
			Light yellow to yellow			
Clarity (post a	autoclaving/heating	g): C	Clear to very slightly opalescent			
(III)Q.C. Test M						
Cultural characteristics observed after an incubation at 35 - 37°C for 18 - 24 hours.						
MICROORGANISM (ATCC)			GROWTH			
Escherichia coli (25922)			Luxuria			
Staphylococcus aureus (25923)			Luxuria			
Enterobacter aerogenes (13048)			Luxuria	int		
Pseudomonas aeruginosa (27853)			Luxuria	int		
Precautions :	1. For Laboratory Use.					
	2. Follow proper, established laboratory procedures in handling and disposing of					
	infectious materials.					
Limitations :	1. Since the nutritional requirements of organisms vary, some strains may be					
	encountered that fail to grow or grow poorly on this medium.					
Use :	For enumeration of microorganisms from water as per ISO6222:1999.					
Storage :	Dehydrated medium- below 30 ° C Prepared mediums- Between 2 to 8°C.					
Packing :	500 gm. bottle					
Product profile:			on	pH (25°C)	Supplement	Sterilization
			ion (500g)		Supplement	Stermzation
B1455	24 g/l	20.83L		7.2 <u>+</u> 0.2	Nil	121ºC/15 min.
				<u> </u>		0, 10
		1		-	1	

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01