

BIOMARK Laboratories – INDIA

www.biomarklabs.com

TECHNICAL SHEET

B237	Mac Conkey Agar w/o CV, NaCl w/ 0.5% Bile Salts	
Formula		
Ingredients :	gms/lit.	
Peptic digest of animal tissue	20.00	
Lactose	10.00	
Bile Salts	5.00	
Neutral red	0.075	
Agar	12.00	
Final pH (at 25°C) : 7.4± 0.2		
Directions :		
Suspend 47.07 grams of medium in 1000 ml distilled water. Heat to boiling with gentle swirling to dissolve the agar completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Avoid overheating. Cool to 45 - 50°C and pour into sterile Petri plates. The surface of the medium should be dry when inoculated. Note: For the cultivation of vibrio species, add 5 grams per litre of sodium chloride before sterilization.		
Principle: Peptic digest of animal tissue is sources of nitrogen and other nutrients. Lactose is a fermentable carbohydrate. When lactose is fermented, a local pH drop around the colony causes a colour change in the pH indicator (neutral red). Bile salts is selective agents that inhibit growth of gram –positive organisms. Agar is a solidifying agent.		
QC Tests – (I)Dehydrated Medium		
Colour :	Light yellow to pink	
Appearance :	Homogeneous Free Flowing powder	
(II)Rehydrated medium		
pH (post autoclaving/heating) :	7.4 ± 0.2	
Colour (post autoclaving/heating) :	Orange red	
Clarity (post autoclaving/heating) :	Clear to slightly opalescent	
(III) Q.C. Test Microbiological		
Cultural characteristics observed after 18 – 24 hours at 35- 37°C.		
MICROORGANISM (ATCC)	GROWTH	COLOUR OF COLONY
Enterobacter aerogenes (13048)	Luxuriant	Pink to red
Escherichia coli (25922)	Luxuriant	Pink to red with bile precipitate
Proteus vulgaris (13315)	Luxuriant	Colourless
Salmonella enteritidis (13076)	Luxuriant	Colourless
Salmonella typhi (6539)	Luxuriant	Colourless
Salmonella paratyphi A	Luxuriant	Colourless
Salmonella paratyphi B	Luxuriant	Colourless
Shigella flexneri (12022)	Fair to good	Colourless
Enterococcus faecalis (29212)	Fair to good	Pale pink to red
Staphylococcus aureus (25923)	Fair to good	Pale pink to red
Precautions :	1. For Laboratory Use.	
	2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.	

Refer disclaimer Overleaf

BIOMARK Laboratories – INDIAwww.biomarklabs.com**TECHNICAL SHEET**

Limitations :	1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
	2. Although MacConkey media are selective primarily for gram – negative enteric bacilli, biochemical and, if indicated, serological testing using pure cultures are recommended for complete identification. Consult appropriate references for further information.				
	3. Incubation of MacConkey Agar plates under increased CO ₂ has been reported to reduce the growth and recovery of a number of strains of gram – negative bacilli.				
Use :	For cultivation and differentiation of enteric bacteria, restricting swarming of Proteus species from specimens like urine which may contain large number of Proteus species as well as potentially pathogenic gram-positive organisms.				
Storage :	Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C.				
Packing :	500 gm. bottle				
Product profile:	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
	B237	47.07g/L	10.622L	7.4± 0.2	NIL

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

BIOMARK Laboratories – INDIA

www.biomarklabs.com

TECHNICAL SHEET