

TECHNICAL SHEET

| | | |
|--|--|--|
| B270I | NITRATE BROTH | |
| Formula | | |
| Ingredients : | gms/lit. | |
| Peptic digest of animal tissue | 5.00 | |
| Meat extract B# | 3.00 | |
| Potassium nitrate | 1.00 | |
| # Equivalent to Beef extract | | |
| Final pH (at 25°C) :7.0 ± 0.2 | | |
| Directions : | | |
| Suspend 9 gms. in 1000 ml. distilled water. Boil to dissolve the medium completely. Dispense in tubes and sterilize by autoclaving at 15 lbs pressure(121°C) for 15 minutes. | | |
| Principles: | | |
| Meat Extract B and Peptone are sources of carbon, protein and nutrients. Potassium Nitrate is a source of nitrate. Nitrate reduction is a valuable criterion for differentiating and identifying various types of bacteria. Certain bacteria reduce nitrates to nitrites only, while others are capable of further reducing nitrite to free nitrogen or ammonia. | | |
| Preparation of Nitrate Test Reagents: 1. Sulfanilic Acid: Dissolve 8 grams of sulfanilic acid in 1 litre 5 N acetic acid. | | |
| 2. Alpha-Naphthylamine reagent: Dissolve 5 grams of alpha-naphthylamine in 1 litre 5 N acetic acid. | | |
| For the test: Add few drops of each reagent i.e. sulphanilic acid (BA078) and a-naphthylamine solution (BA069) into the tube containing culture to be tested. A distinct red or pink colour indicates nitrate reduction. A control (un-inoculated) tube should also be tested. | | |
| QC Tests - (I)Dehydrated Medium | | |
| Colour : | Cream to yellow | |
| Appearance : | Homogeneous Free Flowing powder | |
| (II)Rehydrated medium | | |
| pH (post autoclaving/heating) : | 7.0 ± 0.2 | |
| Colour (post autoclaving/heating) : | Light amber | |
| Clarity (post autoclaving/heating) : | Clear | |
| (III)Q.C. Test Microbiological | | |
| Cultural characteristics observed after 18 - 24 hrs at 35 - 37°C. | | |
| MICROORGANISM (ATCC) | GROWTH | NITRATE REDUCTION |
| Acinetobacter calcoaceticus (19606) | Luxuriant | Negative reaction. |
| Enterobacter aerogenes (13048) | Luxuriant | Positive reaction,distinct red-pinkcolour developed within 1-2 minutes. |
| Escherichia coli (25922) | Luxuriant | Positive reaction,distinct red-pink colour developed within 1-2 minutes. |
| Salmonella typhimurium (14028) | Luxuriant | Positive reaction,distinct red-pink colour developed within 1-2 minutes. |
| Precautions : | <ol style="list-style-type: none"> 1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials. 3. IRRITANT. Irritating to eyes, respiratory system and skin. Avoid contact with skin and eyes. Do not breathe dust. Wear suitable protective clothing. Keep container tightly closed. Target organ(s) : Blood, Nerves. | |
| Limitations : | <ol style="list-style-type: none"> 1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium. 2. The addition of too much zinc dust may cause a false - negative reaction or a momentary colour reaction. 3. The nitrate test is very sensitive. An uninoculated nitrate control should be tested with reagents to determine whether the medium is nitrate free and that the glassware and reagents have not been contaminated with nitrous oxide. 4.The inoculum should not be taken from broth suspension of the organisms. | |
| Use : | For detection of nitrate reduction by bacteria. It is also recommended by ISO committee as per specification ISO 7932:1993 | |

Refer disclaimer Overleaf

BIOMARK Laboratories-INDIAwww.biomarklabs.com**TECHNICAL SHEET**

| | | | | | |
|-------------------------|--|--------------------------------|-----------|------------|--------------------|
| Storage : | Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C. | | | | |
| Packing : | 500 gm. bottle | | | | |
| Product profile: | Reconstitution | Quantity on Preparation (500g) | pH (25°C) | Supplement | Sterilization |
| B270 I | 9g/l | 55.55L | 7.0 ± 0.2 | NIL | 121°C / 15 minutes |

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.