

TECHNICAL SHEET

B347	UREA AGAR BASE (CHRISTENSEN) (AUTOCLAVABLE)		
Formula			
Ingredients :		gms/lit.	
Peptic digest of animal tissue	1.00		
Dextrose	1.00		
Sodium chloride	5.00		
Disodium phosphate	1.20		
Monopotassium phosphate	0.80		
Phenol red	0.012		
Agar	15.00		
Final pH (at 25°C) : 6.8 ± 0.2			
Directions :			
Suspend 24.01 grams in 950 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 10 lbs pressure (115°C) for 20 minutes. Cool to 45-50°C and aseptically add 50 ml of sterile 40% Urea Solution (BF048) and mix well. Dispense into sterile tubes and allow to set in the slanting position. Do not overheat or reheat the medium as urea decomposes very easily			
Principle :			
Peptic digest of animal tissue provides carbon and nitrogen required for good growth of a wide variety of organisms. Dextrose is included as an energy source. Sodium Chloride maintains the osmotic balance of the medium. Potassium phosphate, Monobasic and Sodium Phosphate, Dibasic provide buffering capability. Urea provides a source of nitrogen for those organisms producing urease. This is indicated by a colour change of the pH indicator, Phenol red, from yellow (pH 6.8) to red to pink – red (pH 8.1). Agar is the solidifying agent.			
QC Tests – (I) Dehydrated Medium			
Colour :	Light yellow to light pink		
Appearance :	Homogeneous Free Flowing powder		
(II) Rehydrated medium			
pH (post autoclaving/heating) :	6.8 ± 0.2		
Colour (post autoclaving/heating) :	Yellowish orange		
Clarity (post autoclaving/heating) :	Clear		
(III) Q.C. Test Microbiological			
Cultural characteristics observed on addition of sterile 40% Urea Solution (BF048) after an incubation at 35-37°C for 18-24 hours.			
MICROORGANISM (ATCC)	GROWTH	UREASE	
Enterobacter aerogenes (13048)	Luxuriant	negative reaction, no change	
Escherichia coli (25922)	Luxuriant	negative reaction, no change	
Klebsiella pneumoniae (13883)	Luxuriant	positive reaction, cerise colour	
Proteus vulgaris (13315)	Luxuriant	positive reaction, cerise colour	
Proteus mirabilis (25933)	Luxuriant	positive reaction, cerise colour	
Salmonella typhimurium (14028)	Luxuriant	negative reaction, no change	

Refer disclaimer Overleaf

TECHNICAL SHEET

Precautions :	1. For Laboratory Use.				
	2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :	1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
	2. The alkaline reaction produced in this medium after prolonged incubation may not be caused by urease activity. False positive reactions may occur due to the utilization of peptones (especially in slant agar by <i>Pseudomonas aeruginosa</i> . For example) or other proteins which raise the pH due to protein hydrolysis and the release of excessive amino acid residues. To eliminate possible protein hydrolysis, perform a control test with the same test medium without urea.				
	3. Do not heat or reheat the medium because urea decomposes very easily.				
	4. Urea Agar detects rapid urease activity of only the urease – positive <i>Proteus</i> species. For results to be valid for the detection of <i>Proteus</i> , the results must be read within the first 2 to 6 hours after incubation. Urease – positive <i>Enterobacter</i> , <i>Citrobacter</i> or <i>Klebsiella</i> , in contrast, hydrolyze urea much more slowly, showing only slight penetration of the alkaline reaction into the butt of the medium in 6 hours and requiring 3 to 5 days to change the reaction of the entire butt.				
Use :	Urea Agar Base with the addition of Urea is recommended for the detection of urease production, particularly by members of the genus <i>Proteus</i> .				
Storage :	Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C.				
Packing :	500 gm. bottle				
Product profile:	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
	B347	24.01 g/l	20.824L	6.8 ± 0.2	40% Urea Solution (BF048)

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.