BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BA108 SODIUM HYPOCHLORITE SOLUTION (AVAILABLE CHLORINE 3.5%)

Formula

Preparation:

Sodium hypochlorite 4-6% Water94.21%

Precautions:

- 1. For Laboratory Use.
- 2. Hazardous chemical. Handle carefully.

Use:General disinfection of lab apparatus.

Quality Control:

Appearance: Colourless to pale yellow clear liquid having an odour of chlorine.

Storage: Below 25°C

Packing: 500 ml bottle

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01

Rev: December 2020