

TECHNICAL SHEET

BF007 EGG YOLK TELLURITE EMULSION**FORMULA****Ingredients:**

(100 ml per vial)

Egg Yolk	30.00 ml
Sterile saline	64.00 ml
Sterile 3.5% potassium tellurite solution	6.00 ml
Final pH (at 25°C)	7.6±0.2

Directions:

Warm up the refrigerated Egg Yolk Tellurite Emulsion to 40-45°C. Shake well to attain uniform emulsion (since on refrigeration emulsion has a tendency to form layers or small lumps). Aseptically add 50 ml in 950 ml of sterile, molten Baird Parker Agar Base (B013 & B1432)/Baird Parker Agar w/Sulpha (B940). Mix well and pour into sterile petri plates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: Sterile stabilized tellurite emulsion of egg yolk recommended for identification of Staphylococcus species.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.