BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF047 CC Supplement

FORMULA

Ingredients:

(Per vial sufficient for 500 ml medium)

Cycloheximide 200 mg Chloramphenicol 25 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of distilled water. Mix well and aseptically add it to 500 ml sterile, molten, cooled (45-50°C), Sabouraud Dextrose Agar Base, Modified Emmons (B1110). Mix well and pour into sterile petri plates.

Precautions: 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

Use: Recommended for the selective cultivation of pathogenic fungi.

Storage: Between 2 to 8°C.

Packing: 5 vials