

BF089

PERFRINGENS SUPPLEMENT II

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Oleandomycin phosphate	0.250mg
Polymyxin B sulphate	5000IU

Directions :

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Perfringens Agar Base O.P.S.P. (B671) along with sterile rehydrated contents of 1 vial of Perfringens Supplement I (BF088). Mix well and pour into sterile petriplates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement, recommended for the selective isolation of *Clostridium perfringens*

Storage: Between 2 to 8°C.

Packing : 5 vial