

BF091	T.S.C. SUPPLEMENT (PERFRINGENS T.S.C. SUPPLEMENT)
FORMULA	
Ingredients: (per vial, sufficient for 500 ml medium)	
D-Cycloserine	200.00 mg
Directions:	
Rehydrated the content of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 475 ml of sterile, molten Perfringens agar Base (B282), along with 25 ml of Egg Yolk Emulsion or 500 ml of Perfringens Agar Base (B282)/Perfringes Agar Base (B1496)/TSC Agar Base (B1462)/Tryptose Cycloserine Dextrose Agar Base (B061) or 750 ml of Tryptose Cycloserine Azide Agar Base (B062). Mix well and pour into the sterile Petri plates.	
Precautions :	
1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.	
Use:	An antibiotic supplement recommended for the selective isolation of Clostridium perfringens.
Storage:	Between 2 to 8°C.
Packing:	5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.