

BF098	TRIBUTYRIN
FORMULA	
Ingredients : (Per vial, sufficient for 1000 ml medium)	
Tributyryn	10.00 ml
Directions :	
Add content of 1 vial aseptically to 990 ml of Tributyrin Agar Base (B076) and then sterilize by autoclaving at 15lbs pressure for 15 minutes. Mix well and pour into sterile Petri plates.	
Precautions :	
1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.	
Use :	The supplement is recommended for detection of lipolytic microorganisms.
Storage:	Between 2 to 8°C.
Packing :	10 vial