

BF099

TRICHOMONAS SELECTIVE SUPPLEMENT II

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Streptomycin	250 mg
Penicillin G	500000 Unit

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and avoid frothing of solution. Aseptically add the contents to 460 ml of sterile, molten, cooled (45-50°C) Trichomonas Agar Base (B337) along with 40 ml of inactivated horse serum. Mix well and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antimicrobial supplement recommended for selective isolation of Trichomonas species.

Storage: Between 2 to 8°C.

Packing : 5 vial