BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF123 MYCOPLASMA UROGENITAL SELECTIVE SUPPLEMENT

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Penicillin 100000 Unit

Nalidixic acid 05mg

Amphotericin B 01 mg

Directions:

Aseptically rehydrate the contents of 1 vial in 10 ml sterile distilled water. Aseptically add 1 vial to 500 ml of sterile, cooled (45-50°C) Urogenital Mycoplasma Broth Base (Mycoplasma Urogenital Broth Base)(B1233) along with 1 vial of Vitamino Growth Supplement (BF037) , 1 vial of Urea Solution (BF124) and 50 ml Horse Serum (BF021).Mix well and dispense as desired

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use:An antibiotic supplement recommended for the selective isolation of mycoplasma associated with urogenital infections.

Storage: Store below2-8°C.

Packing: 5 vial

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01

Rev: December 2020