

TECHNICAL SHEET

BF142

MUG SUPPLEMENT (50 MG PER VIAL)

FORMULA

Ingredient :

4-Methylumbiliferyl β -D-Glucuronide (MUG) 50.00 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml methanol. One vial is sufficient for 500 ml agar media or 1000ml of broth media.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.
3. Thawing of the serum should be done as quickly as possible.

Use: A fluorogenic substrate is recommended for measuring β -glucosidase activity for identification of E.coli

Storage: Store between 2-8°C

Packing: vial

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.