

BI002 Agar Powder (Bacteriological)

Specifications

Appearance: Pale cream to cream coloured homogeneous free flowing powder

Clarity (1.5% solution): Clear to slight opalescent gel

pH (1.5% solution) : 6.5 – 7.5

Moisture :NMT 9.0%

Gelling Temperature: 35-38⁰c

Microbial Load:(CFU/g)

Total aerobic microbial count by plate method when incubated at 30-35⁰c for 48 hours

Bacterial count : NMT 5000

Total Yeast and mold count by plate method when incubated at 20-25⁰c for 120 hours

Yeast and mold count :NMT 250

Test for specific organisms

E. coli : Absent

Salmonella species : Absent

Pseudomonas aeruginosa : Absent

Staphylococcus aureus : Absent

Applications:

As solidifying agent in Microbiological culture media at concentrations:

For routine media – 1.2.to 1.7%

For softer media –0.5%

For semisolid media – 0.15%

For inhibition of Proteus spp. – 2.0 %

Precaution

1. For Laboratory Use

2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Storage : Below 30⁰c in cool, dry place.

Packing :500g

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.