

BF060

V. C. N. T. SUPPLEMENT

FORMULA

Ingredients :

(per vial sufficient for 500 ml. of medium)

Vancomycin	1.50 mg
Colistin methane sulphonate	3.75 mg
Trimethoprim	2.50 mg
Nystatin	6,250 units

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml sterile, molten Thayer Martin medium base (B333) / GC Agar Base (B1049) / Transgrow medium base (B078). Mix well and pour into sterile Petri plates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement recommended for the selective isolation of *Nisseria gonorrhoeae* and *Neisseria meningitidis*.

Storage: Between 2 to 8°C.

Packing : 2 vials