

BF146

VANCOMYCIN SUPPLEMENT

FORMULA

Ingredients:

(Per vial, sufficient for 1000 ml medium)

Vancomycin	10 mg
------------	-------

Directions:

Add 10 ml of sterile distilled water in 1 vial. Mix well and aseptically add the contents of 1 vial to 1000ml of sterile, cooled(45-50°C) Modified Lauryl Sulphate Tryptose Broth(B1322) /Modified Lauryl Sulphate Tryptose Broth (B1410). Vancomycin solution can be kept at 0 to 5°C for 15 days)

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective enrichment of *Enterobacter sakazakii*.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.